Prevalence of genomic alterations in Xerna tumor microenvironment subtypes in triple negative breast cancer patients

Gargi D. Basu¹, Janine LoBello¹, Snehal G. Thakkar¹, Jessica Aldrich¹, Nishitha A. Flannery¹, Nishitha Therala¹, David W. Hall¹, Daniel Pointing³, Lea Vohar³, Roman Lustrik³, Luka Ausec³, Mark Uhlik², Seema Iyer², Laura Benjamin², Frederick L. Baehner¹ Exact Sciences¹, OncXerna Therapeutics², Genialis Inc.³

Introduction

Immune checkpoint inhibitors (ICI) are an important therapeutic option for patients with triple negative breast cancer (TNBC). However, identification of patients most likely to respond is challenging. PD-L1 positivity by immunohistochemistry is the standard biomarker used for ICI therapy selection in TNBC. However, other biomarkers, such as analysis of the tumor microenvironment (TME) may be more accurate in predicting response.

The XernaTM TME Panel uses RNA sequencing data and machine learning to analyze the TME, utilizing the angiogenic and immunogenic biology of the TME to classify tumors into four TME subtypes (*Figure* 1). In this study, the distribution of Xerna TME subtypes and associated genomic alterations in TNBC were investigated for their potential use in therapy selection.

Methods

A total of 203 TNBC patient samples underwent tumor-normal whole-exome, whole-transcriptome sequencing testing with the OncoExTraTM assay.

Expression data from whole-transcriptome sequencing were analyzed with the Xerna TME Panel to assign each sample to one of four subtypes:

Immune Active (IA),

Immune Suppressed (IS),


Immune Desert (ID), and

Angiogenic (A).

Actionable alterations, defined as those with FDA-approved matched therapies in any cancer, with matched clinical trials, or with evidence in cancer guidelines or the literature for possible matched therapies, were also identified and associations across Xerna subtypes were explored.

V Figure 1.

The machine learning-based Xerna score is obtained from RNA gene expression levels of ~100 genes. The score reflects the dominant cellular micro-environment of the tumor, along immune and angiogenic axes, and may be useful for predicting response to particular therapies, thus informing therapy decisions. [1]



Subtypes: ID + IA

Suppressed T cell

Inactive T cell

Active T cell

🌲 M1 Macrophage

🕘 M2 Macrophage

PMN

Results

Approximately half (100 of 203; 49.3%) of the patient samples had high (IA+IS) immune subtypes (Table 1). Targetable alterations associated with an FDA-approved therapy were present in 114 (56.2%) patients.

No biomarkers were significantly associated (p<0.05) with high (IA+IS) versus low (ID+A) immune scores (*Table 2*).

Biomarkers associated with ICI response, namely mismatch repair gene alterations (MSH2/3/6, MLH1/3, PMS1/2), high tumor mutational burden (TMB-high) and microsatellite instability were detected in only 6 (3.0%), 3 (1.5%) and 1 (0.5%) patient samples respectively, and all but 1, an MSH6 alteration, were in high immune subtype samples (*Table 2*).

Table 1.

Patient characteristics and Xerna tumor microenvironment subtype / immune group.

			Xerna su	Immune group								
Variable	All Samples	Α	IA	ID	IS	High (IA/IS)	Low (A/ID)					
N	203	48	51	55	49	100	103					
Age (years)												
Mean (SD)	56.6 (13.30)	60.6 (10.68)	53.1 (13.00)	54.7 (13.24)	58.6 (14.90)	55.8 (14.16)	57.5 (12.41)					
Sex												
Female	202 (99.5%)	47 (97.9%)	51 (100.0%)	55 (100.0%)	49 (100.0%)	100 (100.0%)	102 (99.0%)					
Male	1 (0.5%)	1 (2.1%)	(%)	(%)	(%)	(%)	1 (1.0%)					
Actionable Alterations per Sample												
Mean (SD)	2.6 (2.03)	2.4 (1.62)	2.7 (2.48)	3.4 (2.28)	2.0 (1.17)	2.3 (1.97)	2.9 (2.05)					
Median	2.0	2.0	2.0	3.0	2.0	2.0	2.0					
Q1-Q3	1 - 3	1 - 3	1 - 3	2 - 5	1 - 3	1 - 3	2 - 4					
Min, Max	0, 16	0, 7	1, 16	1, 11	1, 5	1, 16	0, 11					

Questions? Please contact gbasu@exactsciences.com

This presentation is the intellectual property of the author/presenter. Contact them for permission to reprint and/or distribute.

Table 2.

Frequency of the 6 actionable biomarkers that were present at least 10 TNBC samples as well as 6 biomarkers associated with ICI response. No correction for multiple comparisons was employed. N is the number of patient samples.

	Xerna subtype				Immune group						
Total (n=203)	A (n=48)	IA (n=51)	ID (n=55)	IS (n=49)	High (IA/IS) (n=100)	Low (A/ID) (n=103)	p-value				
173 (85.2%)	36 (75.0%)	48 (94.1%)	49 (89.1%)	40 (81.6%)	88 (88.0%)	85 (82.5%)	0.32				
34 (16.7%)	8 (16.7%)	9 (17.6%)	9 (16.4%)	8 (16.3%)	17 (17.0%)	17 (16.5%)	1.00				
17 (8.4%)	5 (10.4%)	1 (2.0%)	5 (9.1%)	6 (12.2%)	7 (7.0%)	10 (9.7%)	0.61				
16 (7.9%)	3 (6.3%)	4 (7.8%)	6 (10.9%)	3 (6.1%)	7 (7.0%)	9 (8.7%)	0.80				
12 (5.9%)	0 (0.0%)	4 (7.8%)	7 (12.7%)	1 (2.0%)	5 (5.0%)	7 (6.8%)	0.77				
12 (5.9%)	5 (10.4%)	3 (5.9%)	1 (1.8%)	3 (6.1%)	6 (6.0%)	6 (5.8%)	1.00				
Biomarkers relevant to ICI therapy											
5 (2.5%)	0 (0.0%)	2 (3.9%)	1 (1.8%)	2 (4.1%)	4 (4.0%)	1 (1.0%)	0.21				
1 (0.5%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	1 (2.0%)	1 (1.0%)	0 (0.0%)	0.49				
1 (0.5%)	0 (0.0%)	1 (2.0%)	0 (0.0%)	0 (0.0%)	1 (1.0%)	0 (0.0%)	0.49				
3 (1.5%)	0 (0.0%)	2 (3.9%)	0 (0.0%)	1 (2.0%)	3 (3.0%)	0 (0.0%)	0.12				
1 (0.5%)	0 (0.0%)	1 (2.0%)	0 (0.0%)	0 (0.0%)	1 (1.0%)	0 (0.0%)	0.49				
	(n=203) 173 (85.2%) 34 (16.7%) 17 (8.4%) 16 (7.9%) 12 (5.9%) 12 (5.9%) relevant to ICI 5 (2.5%) 1 (0.5%) 1 (0.5%) 3 (1.5%)	(n=203) $(n=48)$ $173 (85.2%)$ $36 (75.0%)$ $34 (16.7%)$ $8 (16.7%)$ $17 (8.4%)$ $5 (10.4%)$ $16 (7.9%)$ $3 (6.3%)$ $12 (5.9%)$ $0 (0.0%)$ $12 (5.9%)$ $5 (10.4%)$ $elevant to ICI therapy$ $5 (2.5%)$ $0 (0.0%)$ $1 (0.5%)$ $0 (0.0%)$ $1 (0.5%)$ $0 (0.0%)$ $3 (1.5%)$ $0 (0.0%)$	Total (n=203)A (n=48)IA (n=51)173 (85.2%)36 (75.0%)48 (94.1%)34 (16.7%)8 (16.7%)9 (17.6%)17 (8.4%)5 (10.4%)1 (2.0%)16 (7.9%)3 (6.3%)4 (7.8%)12 (5.9%)0 (0.0%)4 (7.8%)12 (5.9%)5 (10.4%)3 (5.9%)12 (5.9%)0 (0.0%)2 (3.9%)elevant to ICI Herapy5 (2.5%)0 (0.0%)1 (0.5%)0 (0.0%)1 (2.0%)3 (1.5%)0 (0.0%)2 (3.9%)	Total (n=203)A (n=48)IA (n=51)ID (n=55)173 (85.2%) $36 (75.0\%)$ $48 (94.1\%)$ $49 (89.1\%)$ $34 (16.7\%)$ $8 (16.7\%)$ $9 (17.6\%)$ $9 (16.4\%)$ $17 (8.4\%)$ $5 (10.4\%)$ $1 (2.0\%)$ $5 (9.1\%)$ $16 (7.9\%)$ $3 (6.3\%)$ $4 (7.8\%)$ $6 (10.9\%)$ $12 (5.9\%)$ $0 (0.0\%)$ $4 (7.8\%)$ $7 (12.7\%)$ $12 (5.9\%)$ $5 (10.4\%)$ $3 (5.9\%)$ $1 (1.8\%)$ relevant to ICI therapy $5 (2.5\%)$ $0 (0.0\%)$ $2 (3.9\%)$ $1 (1.8\%)$ $1 (0.5\%)$ $0 (0.0\%)$ $1 (2.0\%)$ $0 (0.0\%)$ $3 (1.5\%)$ $0 (0.0\%)$ $2 (3.9\%)$ $0 (0.0\%)$	Total (n=203)A (n=48)IA (n=51)ID (n=55)IS (n=49)173 (85.2%) 36 (75.0%) 48 (94.1%) 49 (89.1%) 40 (81.6%) 34 (16.7%) 8 (16.7%) 9 (17.6%) 9 (16.4%) 8 (16.3%) 17 (8.4%) 5 (10.4%) 1 (2.0%) 5 (9.1%) 6 (12.2%) 16 (7.9%) 3 (6.3%) 4 (7.8%) 6 (10.9%) 3 (6.1%) 12 (5.9%) 0 (0.0%) 4 (7.8%) 7 (12.7%) 1 (2.0%) 12 (5.9%) 5 (10.4%) 3 (5.9%) 1 (1.8%) 3 (6.1%) 12 (5.9%) 0 (0.0%) 2 (3.9%) 1 (1.8%) 2 (4.1%) 1 (0.5%) 0 (0.0%) 1 (2.0%) 0 (0.0%) 1 (2.0%) 1 (0.5%) 0 (0.0%) 1 (2.0%) 0 (0.0%) 1 (2.0%) 3 (1.5%) 0 (0.0%) 2 (3.9%) 0 (0.0%) 1 (2.0%)	Total (n=203)A (n=48)IA (n=51)ID (n=55)IS (n=49)High (IA/IS) (n=100)173 (85.2%)36 (75.0%)48 (94.1%)49 (89.1%)40 (81.6%)88 (88.0%)34 (16.7%)8 (16.7%)9 (17.6%)9 (16.4%)8 (16.3%)17 (17.0%)17 (8.4%)5 (10.4%)1 (2.0%)5 (9.1%)6 (12.2%)7 (7.0%)16 (7.9%)3 (6.3%)4 (7.8%)6 (10.9%)3 (6.1%)7 (7.0%)12 (5.9%)0 (0.0%)4 (7.8%)7 (12.7%)1 (2.0%)5 (5.0%)12 (5.9%)5 (10.4%)3 (5.9%)1 (1.8%)3 (6.1%)6 (6.0%)12 (5.9%)5 (10.4%)2 (3.9%)1 (1.8%)2 (4.1%)4 (4.0%)1 (0.5%)0 (0.0%)0 (0.0%)0 (0.0%)1 (2.0%)1 (1.0%)1 (0.5%)0 (0.0%)1 (2.0%)0 (0.0%)1 (1.0%)3 (3.0%)3 (1.5%)0 (0.0%)2 (3.9%)0 (0.0%)1 (2.0%)3 (3.0%)	Total (n=203)A (n=48)IA (n=51)ID (n=55)IS (n=49)High (IA/IS) (n=100)Low (A/ID) (n=103)173 (85.2%)36 (75.0%)48 (94.1%)49 (89.1%)40 (81.6%)88 (88.0%)85 (82.5%)34 (16.7%)8 (16.7%)9 (17.6%)9 (16.4%)8 (16.3%)17 (17.0%)17 (16.5%)17 (8.4%)5 (10.4%)1 (2.0%)5 (9.1%)6 (12.2%)7 (7.0%)10 (9.7%)16 (7.9%)3 (6.3%)4 (7.8%)6 (10.9%)3 (6.1%)7 (7.0%)9 (8.7%)12 (5.9%)0 (0.0%)4 (7.8%)7 (12.7%)1 (2.0%)5 (5.0%)7 (6.8%)12 (5.9%)5 (10.4%)3 (5.9%)1 (1.8%)3 (6.1%)6 (6.0%)6 (5.8%)12 (5.9%)5 (10.4%)2 (3.9%)1 (1.8%)2 (4.1%)4 (4.0%)1 (1.0%)10 (0.5%)0 (0.0%)0 (0.0%)1 (2.0%)1 (1.0%)0 (0.0%)1 (0.5%)0 (0.0%)1 (2.0%)0 (0.0%)1 (1.0%)0 (0.0%)1 (0.5%)0 (0.0%)1 (2.0%)0 (0.0%)1 (1.0%)0 (0.0%)3 (1.5%)0 (0.0%)2 (3.9%)0 (0.0%)1 (2.0%)3 (3.0%)0 (0.0%)				

Conclusions

- to ICI therapy.
- considered for novel combination therapies.

Reference

1. Uhlik M, Pointing D, Iyer S, Ausec L, Štajdohar M, Cvitkovič R, Žganec M, Culm K, Santos VC, Pytowski B, Malafa M, Liu H, Krieg AM, Lee J, Rosengarten R, Benjamin L. Xerna[™] TME Panel is a machine learning-based transcriptomic biomarker designed to predict therapeutic response in multiple cancers. Front Oncol. 2023 May 12;13:1158345. doi: 10.3389/fonc.2023.1158345. PMID: 37251949; PMCID: PMC10213262.

Acknowledgements

The study was funded by Exact Sciences. Special thanks to Karen White for assistance in poster preparation.



• The Xerna TME Panel classified 49.3% of TNBC patient tumors to IA or IS, suggesting they may respond

• Many (56.2%) patient tumors harbored alterations associated with potential targeted therapies (e.g., PIK3CA and PTEN with PI3K/mTOR inhibitors and BRCA1 with PARP inhibitors), which might be

• These findings warrant further study and clinical validation in TNBC patients treated with ICI therapy.